


Wireless Scale

Model # ADF-B833



Contents

INTRODUCTION.....	3
NOTES ON SAFETY	3
⚠PRECAUTION NOTICE.....	3
⚠IMPORTANT NOTICE TO USERS.....	3
⚠CARE AND MAINTENANCE	4
PRODUCT DESCRIPTION.....	5
Front View	5
LED Readout.....	5
Back View	6
Measurement Unit Switch.....	6
PREPARATION BEFORE USE	7
Insert the Batteries.....	7
Remove the Cling Label.....	8
Set the Scale On The Floor	8
Start the Scale	8
OPERATION	9
Foot Placement	9
WEIGHING YOURSELF.....	10
TO TURN OFF THE SCALE	10
TROUBLESHOOTING	11
CARE AND MAINTENANCE.....	11
SPECIFICATIONS	13
CONTACT INFORMATION AND NOTICES.....	14
FCC Compliance Statement	15

Wireless Scale
INSTRUCTION MANUAL Model # ADF-B883

INTRODUCTION

This scale is provided to you by Remote Care Partners (RCP), for use exclusively with the health monitoring program from your provider.

Please read this manual carefully and keep it handy for future reference.

NOTES ON SAFETY

Please read this section carefully to familiarize yourself with features and operations before using the unit.

- The warning signs and the sample icons shown here are listed for you to use this product safely and correctly as well as to prevent product damage, risk, and injury to you and others.
- The icons and meanings are as follows:

 PRECAUTION NOTICE	Indicates the right conditions to use the product and to prevent damage, risk, and injury.
 IMPORTANT NOTICE TO USERS	Indicates important notices users should read before using the product.

CARE AND MAINTENANCE

Indicates matters in which the possibility of damage may happen as a result of incorrect handling and improper maintenance.

Precaution Notice


OVERLOAD WARNING: the maximum weight capacity of this scale is 330 lbs/150 kgs. Remove the weight immediately when the scale displays “Err”; otherwise, permanent damage to the scale will occur.

This product is not for infants and small children.


- Make sure to use only the type of battery stated (see Section “PREPARATION BEFORE USE”).
- For accurate results, place your feet in the same position every time you weigh yourself.
- Stand still during the measurement
- Weigh yourself without clothing whenever possible. Clothes can vary in weight.
- For the most consistent reading, weigh yourself at the same time of day.

PRODUCT DESCRIPTION


Front View



LED Readout



Back View



Measurement Unit Switch


You can change the measuring unit between kg and lb by pressing the button located on the bottom of the scale.

PREPARATION BEFORE USE

Insert the Batteries

New batteries were installed in your device before shipping. Unless instructed by a support person, skip this section

This scale operates on 3 AAA batteries (included). Remove any plastic wrap from the batteries before installation. Always replace both batteries at the same time; do not mix old and new batteries. Open the battery cover on the back of the scale unit. Place the batteries into the battery compartment according to the "+" and "-" polarity symbols at the bottom of the compartment. Replace the battery compartment cover. If you do not intend to use this unit for a prolonged period of time, it is advisable to remove the batteries before storing.



TIPS

**Replace batteries when
"Lo" is displayed**

Remove the Cling Label

Some models have a static cling label on the LED lens to prevent scratching. Remove before use.

Set the Scale On The Floor

Always place the scale on a flat, hard surface. Do not use on carpets

Start the Scale


Press lightly on the scale platform to start. The Scale Display will show 0.0 and then the display will turn off. You can now step on the scale to get a measurement

OPERATION


The accuracy of the results depends on how you stand on the scale. Position your feet square to the scale, parallel and equally distant from the outer edges of the scale.

Foot Placement

It is important to get your feet placed correctly on the scale. Incorrect foot placement will result in an inaccurate measurement



1 Correct Foot Placement




2 Incorrect Foot Placement

Do not allow other people to use your scale. The scale is specifically assigned to you and you only. Allowing other people to use your scale can mislead your provider.

WEIGHING YOURSELF

1. Step on the scale (disregard the reading) and then step off. This process turns on the scale and makes it ready for use.
2. Step onto the scale and stand still while the weight is being measured.
3. The screen displays the weight.
4. The result will automatically be transferred to the RCP Hub and sent to your provider.
5. The scale turns off automatically after use.



TO TURN OFF THE SCALE

The scale turns off automatically after a few seconds if no buttons are pressed or nothing is weighed.

TROUBLESHOOTING

No weight displayed

Check to see if the scale is powered up and started up. If not, please refer to the Section "Preparation Before Use".

Why does the display read "Lo"?

Battery is running low. Replace the batteries.

CARE AND MAINTENANCE

- **OVERLOAD WARNING:** the maximum weight capacity of this scale is 330 lbs/150 kgs. Remove the weight immediately when the scale displays "Err"; otherwise, permanent damage to the scale will occur.
- The product is intended for home/consumer use only; it is not intended for professional use in hospitals or medical facilities.
- Do not disassemble the product. Other than replacing the batteries, it contains no user-serviceable parts.
- Clean after use with a lightly dampened cloth. Do not use solvents or immerse the product in water.
- Your scale contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
- Do not store the scale where you keep cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store the scale on its side.
- This scale is a sensitive weighing device. To prevent running down the batteries, do not store

anything on the scale.

- Do not dispose of batteries in a fire. Batteries may explode or leak. **Remove the batteries if the scale will not be used for a long period of time.**
- When replacing batteries, if there is more than one battery in the scale, replace all batteries at the same time; do not mix old and new batteries. Do not mix alkaline, standard (carbon-zinc), or rechargeable (ni-cad, ni-mh, etc) batteries.

NOTE: Please recycle or dispose of batteries per local regulations.

WARNING: Batteries may pose a choking hazard. As with all small items, do not let children handle batteries. If swallowed, seek medical attention immediately.

PRECAUTION: Do not dispose of batteries in a fire. Remove the battery if the scale will not be used for a long period of time.

Caution: Changes or modifications not expressly approved by the partner responsible for compliance could void the user's authority to operate the equipment.

SPECIFICATIONS

Dimensions	Scale 300 x 300 x25MM
	Glass 300 x 300 x6mm
Net Weight	~1.6kg
Display	Digital LCD with Blue Back
LED Size	85 x 45 x 2.8 mm
Measurement Unit	Kilogram/Stone/Pound
Measurement Range	5kg to 150 kg
	0st: 11lb to 23st: 8pb
	11lb to 30lb
Increment	0.1kg / 0.2 lb
Accuracy	5-50Kg $\pm$ 0.4Kg
	100kg $\pm$ 0.6Kg
	150kg $\pm$ 0.8Kg
Working Environment	Temp: 5°C-40°C 40°F-104°F Humidity: <90%
Storage Environment	Temp: -20°C-60°C -4°F-140°F
Auto-On	SENSE ON technology
Auto-OFF	Automatic

CONTACT INFORMATION AND NOTICES

Remote Care Partners

855-477-7000

www.remotecarepartners.com

AndseFit LTD

Unit 513 Lakeside 1,

No. 8 Science Park West Ave

Hong Kong Science Park

Hong Kong

FCC Compliance Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

This device may not cause harmful interference, and

This device must accept any interference received, including interference that may cause undesired operation.

FCC WARNING

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver.
- Connect the equipment into an outlet different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.